

# TruCLEAN™ Stainless Steel Dosage Sieves

Precisely positioned apertures



30-46/30-47  
TruCLEAN  
Stainless Steel Sieves

## Description

Stainless Steel Dosage Sieves mount conveniently on 36-Liter bucket or waste containment bucket, depending on sieve size. It's precisely positioned apertures release excess disinfectant into the waste containment bucket. TruCLEAN Deluxe Disinfection System comes equipped with one Stainless Steel Sieve (No. 30-46) for use with Waste Containment Bucket (No. 30-09) and one Sieve (No. 30-47) for use with 36-Liter Bucket (No. 30-36). Both options are designed for small area cleaning and disinfecting. Compatible with gamma, ETO autoclave sterilization at 121° C for 30 minutes.

#30-46 Dimensions: 16 3/4" length x 5 3/8" width x 1 1/8" height

#30-47 Dimensions: 17" length x 5 1/8" width x 1 3/8" height

## Features & Benefits

- Optimize application of disinfecting solutions
- New design allows for use with waste containment bucket
- Autoclavable, gamma ETO sterilization compatible
- Impervious to most chemicals, will not rust, corrode or oxidize

## Applications

- Compatible with TruCLEAN Deluxe Disinfection Systems
- Designed for small area cleaning and disinfecting

## Product Information

Part #	Description	Material
30-46	TruCLEAN Stainless Steel Sieve (For use with (30-09) waste containment)	Type 304 Stainless Steel
30-47	TruCLEAN Stainless Steel Sieve (For use with (30-36) 36-Liter bucket)	Type 304 Stainless Steel

## Ordering Information

Part #	Description	Carton Size (L x W x H)	Carton Weight	Packaging
30-46	TruCLEAN Stainless Steel Sieve	18" x 13" x 4"	3 lbs.	1/each
30-47	TruCLEAN Stainless Steel Sieve	18" x 13" x 4"	3 lbs.	1/each

Notes: Data shown reflects typical values applicable to the part numbers listed. Please contact a Perfex sales representative for more details.



32 Case Street  
Poland, NY 13431  
Tel: 800-848-8483  
315-826-3600  
Fax: 315-826-7471  
[www.PerfexOnline.com](http://www.PerfexOnline.com)  
e-mail: [perfex@perfexonline.com](mailto:perfex@perfexonline.com)

