

DESCRIPTION AND GENERAL PROPERTIES

- Hand specific design provides optimum comfort, reduced hand fatigue, available in sizes 6.5 - 9 including half sizes.
- Good tactile sensitivity and good resistance

MAIN FIELDS OF USE

Mapa Advantech

- Controlled Environment: Product protection

Electronics

- Photomask assembly, integrated circuit assembly, CD production, placement of thin layers

Health

- Pharmaceutical manufacturing, manufacturing of medical equipment

Laboratory

- Research, analysis, manipulation of precision parts, optical processes

INSTRUCTIONS FOR USE

Instructions for use

- It is not recommended to persons sensitized to natural latex, dithiocarbamates and thiazoles to use these gloves .
- Put the gloves on dry, clean hands.
- Position the cuff over the garment to prevent skin exposure.
- These gloves are designed for single use only.
- If gloves have been contaminated with chemicals, do not touch the external side of the gloves when taking them off. Fold back the cuff end and pull them off while turning them inside out.

Storage conditions

Store the gloves in their sealed packaging protected from direct sunlight, far from heat sources and electrical equipment

Food contact US

FDA 21CFR 177.2600

LEGISLATION

This product is not classified as hazardous according to the directive 1999/45/EC of the European Parliament and of the Council. This product does not contain any substances of very high concern according to the regulation n° 1907/2006 of the European Parliament and of the Council (REACH).

- CE-Type Examination Certificate

0075/014/162/05/05/0074

- Issued by the notified body nr

0075

C.T.C – F-69367 LYON CEDEX 07

- Quality assurance certificat

0334 Asqual 14 rue des Reculettes -F-75013 PARIS