



KIMTECH PURE* M3 Sterile Pleat-Style Face Mask

Codes: 62467, 62470, 62494

Product Information

Code	Description	Color	Style	Size	Case Count
62467	M3 Sterile Face Mask with Soft Ties	White	Pleated	7"	20/bag, 200/case
62470	M3 Sterile Face Mask with Knitted Earloops	White	Pleated	7"	20/bag, 200/case
62494	M3 Sterile Face Mask with Soft Ties	White	Pleated	9"	20/bag, 200/case



Physical Properties	Results	Test Method
Particle Filtration Efficiency @ 0.1 micron, Post-Sterilization	89.3%	ASTM F1215-89
Differential Pressure in mm H ₂ O @ 8 LPM flow rate, Post-Sterilization	2.16	MIL-M-36954C
Bacterial Filtration Efficiency @ 3.0 micron, Post-Sterilization	94.1%	MIL-M-36954C



Mask Construction	
Binding / Ties	White, Spunlaced Polyester
Nosepiece	Fully enclosed, soft, malleable aluminum
Outer Facing	Clear, Polyethylene, Apertured Film
Inner Facing	White BiCoSof* Fabric ¹
Filter Media	White, Meltblown Polyester
Construction	Ultra Sonic Bonding

¹BiCoSof Fabric is a Spunbond Material made from a bicomponent polypropylene/polyethylene (PP/PE) fiber.

Additional Information:

- Manufactured in ISO 9001:2000 registered facilities.
- All Face Masks are double bagged.
- Recommended for Class 3 or higher cleanroom environments.
- Gamma irradiated with Sterility Assurance Level of 10⁻⁶.
- Certificates of Irradiation are available online at www.kimtech.com/certificates.

Workmanship Requirements

Face Mask shall be made in accordance with reasonable industry practice with respect to holes, tears, wrinkles, cleanliness, odors and foreign materials or dirt.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers and are warranted for the following periods from end user's date of purchase (verified by valid sales receipt) (a) five years for Balder Technology auto darkening filters; (b) two years for all other auto darkening filters; and (c) one year for powered air purifying respirators; (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

©*Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. © 2014 KCWW.

