



G5 Sterile Latex[†] 12" Gloves Hand Specific

Codes: 56851, 56852, 56853, 56854, 56855, 56856, 56857, 56850

Product Information

Code	Affiliate #	Size	Palm Width*(mm)	Length* (mm)	Quantity
56851	HC1160S	6.0	77	290	20 pairs/bag; 10 bags/case
56852	HC1165S	6.5	83	305	20 pairs/bag; 10 bags/case
56853	HC1170S	7.0	89	305	20 pairs/bag; 10 bags/case
56854	HC1175S	7.5	95	305	20 pairs/bag; 10 bags/case
56855	HC1180S	8.0	102	305	20 pairs/bag; 10 bags/case
56856	HC1185S	8.5	108	305	20 pairs/bag; 10 bags/case
56857	HC1190S	9.0	114	305	20 pairs/bag; 10 bags/case
56850	HC1110S	10.0	130	305	20 pairs/bag; 10 bags/case



*As measured by ASTM D 3577; 2.5 AQL¹

Property	Description
General Description	Natural rubber latex, powder free, hand specific gloves, textured palm and palm side of fingers, curved finger design, beaded cuff, silicone free, double chlorination

Physical Property	Unit	Upper Limit	Test Method (1.5 AQL)
Watertightness	NA	Pass	ASTM D 5151
Sterility	SAL	10 ⁻⁶	ANSI/AAMI/ISO 11137

Thickness	Nominal	Nominal	Test Method (2.5 AQL)
Middle Finger	8.66 mil	0.22 mm	ASTM D 3577
Palm	7.87 mil	0.20 mm	ASTM D 3577
Cuff	5.51 mil	0.14 mm	ASTM D 3577

Property	Tensile Strength	Ultimate Elongation	Median Force at Break	Protein Content (Max)
Before Aging	24 Mpa	750%	10.5 Newtons	
After Aging	18 Mpa	560%	7.5 Newtons	50 µg/g
Test Method (2.5 AQL)	ASTM D 412 ASTM D 573	ASTM D 412 ASTM D 573	EN 455-2	ASTM D 5712

¹ AQL levels defined according to ISO 2859-1 (current edition) Sampling Procedures for Inspection by Attributes.

† This product contains natural rubber latex, which may cause allergic reactions. Safe use of this glove by or on latex sensitive individuals has not been established. This product should not be worn, or exposed to, individuals allergic to natural rubber latex.

Workmanship Requirement:

Gloves shall be made in accordance with reasonable industry practice with respect to defects, dirt and contamination. Gloves should be packaged as specified having described size and color in the appropriate quantity as listed. Shipping Case and inner packaging shall be labeled with product name, size, lot/batch and catalog number.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers and are warranted for the following periods from end user's date of purchase (verified by valid sales receipt) (a) five years for Balder Technology auto darkening filters; (b) two years for all other auto darkening filters; and (c) one year for powered air purifying respirators; (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

©*/Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. © 2013 KCWW.



Exceptional Workplaces

Rev: Nov 2013