

# Undercounter and Countertop Freezers

Free Standing Freezers with Solid Door

## UCF Series

Door type Solid

Internal capacity 42L / 1.5 cu.ft. 70L / 2.5 cu.ft. 113L / 4 cu.ft.

### Specifications

Model	UCF000	UCF001	UCF101A	UCF101B
Capacity	42L / 1.5 cu.ft.	70L / 2.5 cu.ft.	113L / 4 cu.ft.	113L / 4 cu.ft.
Door	Single solid, locking			
Temperature range	-20°C / -4°F	0°C to -20°C / +32°F ~ -4°F	0°C to -30°C / +32°F to -22°F	-30°C to -40°C / -22°F to -40°F
Temperature control	Mechanical thermostat	Digital display with Hi/Low alarm, Alarm relay dry contacts		
Alarm system	-	Emits an audible and visual signal when the temperature gets to hi or low. The alarm has an over and under temperature setting, alarm silencing switch, and dry contact relay. A relay for a remote alarm hook-up is also provided.		
Access port	Probe access port 3/8"	1/2"		1/2" rear
Shelves		2 full	2 adjustable, 1 full, 1 half	3 adjustable
Insulation	High density urethane foam	Polyurethane foam		
Refrigerant	R134a CFC-Free refrigerant	R600		R290
Compressor	-	-	1/3 hp	1/4 hp
Air circulation	-	Forced air		
Defrost	Manual			
Exterior construction	Painted white			
Internal dimensions (WxDxH)	14" x 13.25" x 14.5" 35 x 33 x 36 mm	17" x 15" x 17" 43 x 38 x 43 cm	19.5" x 15.75" x 22" 50 x 40 x 56 cm	20" x 19.5" x 24.25" 50 x 49 x 61 cm
External dimensions (WxDxH)	18.5" x 19.5" x 19.5" 46 x 49 x 49 cm	20" x 21.5" x 25" 50 x 54 x 63 cm	23.75" x 21.5" x 32" 60 x 64 x 84 cm	24" x 24" x 33" 61 x 61 x 83 cm
Weight	50 lbs.	110 lbs.	150 lbs.	
Voltage	115V, 60 Hz, 1 phase		115V, 60 Hz, 1 phase	
Running amps / breaker	1.1A, 310 BTUs/hr	1.2A / 15A	2 / 15	5 / 15
Supply plug	-		5-15P NEMA	
Optional accessory	Chart recorder			



42L  
UCF001



70L  
UCF001



113L  
UCF101A

113L  
UCF101B