

TECHNICAL DATA SHEET

Part#: 36PH

Hemostats - Curved ★★- SS - Vinyl-Coated handle - 5"

The 36PH two star stainless steel locking hemostat with 25°curved serrated jaws has vinyl-coated handles. The overall length is 5".

**

SS (AISI 440)

Insulative

Curved

2

• Precision Level

Material

Surface Resistivity

Shape

• Group Small Hand Tools

• Angle 25°

Body Width 0.26" (6.5mm)
Body Height 0.172" (4.3mm)
OAL 5" (125mm)

• Weight (OZ)

UPC
Schedule B
ECCN
UNSPSC Code
Country of Origin
662847003851
9018.90.80.00
EAR99
42151685
Pakistan





