

TECHNICAL DATA SHEET

Part#: 38PH

Hemostats - Curved -★★- SS - Vinyl-Coated Handle - 6"

The 38PH two star stainless steel locking hemostat with 30° curved serrated jaws has vinyl-coated handles. The overall length is 6".

• Precision Level	★★
• Material	SS (AISI 440)
• Surface Resistivity	Insulative
• Shape	Curved
• Group	Small Hand Tools
• Angle	30°
• Body Width	0.283" (7.08mm)
• Body Height	0.193" (4.83mm)
• OAL	6" (150mm)
• Weight (OZ)	2
• UPC	662847003912
• Schedule B	9018.90.80.00
• ECCN	EAR99
• UNSPSC Code	42151685
• Country of Origin	Pakistan

