

Visit us online now:
torreypinsscific.com

760-930-9400

Toll Free within the USA:
866-573-9104

PRODUCT DATA SHEET

IC30 AND IC35

ECHOTHERM™ CHILLING/HEATING DRY BATHS

EchoTherm Models IC30 Digital and IC35 Programmable High Capacity Series Chilling/Heating Dry Baths feature temperature ranges from -20°C to 100°C



- Model IC30 (right)
- Model IC30 with attachments (left)

- Peltier powered for chilling and heating
- Settable temperature ranges (useable range will vary with block size):
 - IC30/IC35: -10°C to 100°C
 - IC30XT/IC35XT: -20°C to 100°C
- Peltier chilling/heating power
 - IC30/IC35: 100 watts
 - IC30XT/IC35XT: 100 watts
- 1°C readable, settable and controllable
- High capacity. For example they use blocks that holds 64- 1.5 ml centrifuge tubes or 9- 50 ml centrifuge tubes
- IC35 models fully programmable with 5-program memory
- Backlit two-line alphanumeric display with 16 characters per line

- 30-day count down timer with alarm and user settable Auto-Off
- RS232 I/O port for control or data collection
- Built-in data logger
- More than 30 standard sample blocks for centrifuge tubes, assay plates, vials, test tubes, round-bottom flasks, and more available as accessories
- Custom sample blocks made on special order for most any use
- Bench top power supply provided for use anywhere in the world
- UL, CSA and CE compliant
- 12-month warranty

Some Applications

- Enzyme reactions and deactivations
- Ligation reactions
- Hybridizations
- PCR sample preparation
- Storing oocytes
- Storing enzymes and DNA libraries
- Incubating samples above, below and at room temp.
- Storing samples at 4°C without an ice bucket
- General molecular biology



TORREY PINES®
SCIENTIFIC
Degrees Better

2713 Loker Avenue West, Carlsbad, CA 92010
Tel: 760-930-9400 | Toll Free within the USA: 866-573-9104
Fax: 760-930-9480
info@torreypinsscific.com

Visit us online now:
torreypinsscific.com

760-930-9400

Toll Free within the USA:
866-573-9104

PRODUCT DATA SHEET

IC30 AND IC35

ECHOTHERM™ CHILLING/HEATING DRY BATHS

EchoTherm Models IC30 Digital and IC35 Programmable High Capacity Series Chilling/Heating Dry Baths feature temperature ranges from -20°C to 100°C



DIMENSIONS

8.5" (216 mm) wide x 10.0" (254 mm)
front to back x 4.0" (102 mm) tall

Weight unpacked: 6 pounds (2.7 kgs)

Weight packed: 11 pounds (5 kgs)