

AQUAHOL STERILE ALCOHOL

Product Specification

| | |
|------------------------------|--|
| Product Name: | AquaHol Sterile Alcohol |
| Manufacturer Part Number(s): | IC-220 |
| Description: | <p>Sterile hand and glove sanitizer comprised of 70% isopropyl alcohol with 30% WFI quality purified water, filtered to 0.1 microns and meets USP requirements. Ideal for use over gloves prior to entering into aseptic areas and between critical processes. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.</p> |
| Manufacturing Standards: | <p>ISO 9001 Certified Facility. 10⁻⁶ Sterility Assurance Level per VDMax Sterilization Method</p> |
| Physical Data: | |
| Materials of Construction: | |
| Sanitizer: | 70% Isopropyl alcohol with 30% purified water/WFI quality, filtered to 0.1 microns. |
| Container: | HDPE |
| Dimensions: | |
| Net Contents | 500 mL |
| For Use with: (as needed) | For use with the ICDISP-3 Touch-free MicroDispenser |
| Autoclavability: | Not Autoclavable |
| Chemical Compatibility: | The product is chemically stable. Under normal conditions of storage and use, hazardous reactions will not occur. |
| Shelf Life: | 3 years after manufacture date |
| Packaging: | 6 bottles per poly-lined case |
| Traceability: | Lot numbers printed on each case and on each individual bottle |
| Irradiation Standards: | Irradiated |
| Certificates: | Certificate of Analysis provided with each order and available upon request |