

# PRODUCT SPECIFICATIONS



## GRPR-(SIZE)-1-T8

Double Chlorinated



PACKED 100 GLOVES PER BOX, 10 BOXES PER CASE (1,000 GLOVES PER CASE)

### NATURAL, PREMIUM, POWDER FREE, LATEX

MATERIAL	CUFF THICKNESS	PALM THICKNESS	FINGER THICKNESS	LENGTH	COLOR
Latex: Made with 100% Natural Latex	4.5 mil +/- 0.5 mil	8.0 mil +/- 0.5 mil	10.0 mil +/- 0.5 mil	9.65" +/-0.25"	NATURAL

**ACCEPTANCES:** All Safety Zone glove components comply with the provisions of the Federal Food, Drug, and Cosmetic Act which allow repeated use in direct contact with food

### PACKAGING SPECIFICATIONS

MFG#	SIZE	GLOVE WEIGHT* (grams)	CASE WEIGHT (lbs)	UPC#	Pallet	HIGH	TIE	CUBIC FEET	COUNTRY OF ORIGIN
GRPR-SM-1-T8	SMALL	9.8	26.5	76358330191	40	4	10	1.18	INDONESIA THAILAND MALAYSIA
GRPR-MD-1-T8	MEDIUM	10.0	27.0	76358330190	40	4	10	1.18	INDONESIA THAILAND MALAYSIA
GRPR-LG-1-T8	LARGE	11.0	29.5	76358330189	40	4	10	1.18	INDONESIA THAILAND MALAYSIA
GRPR-XL-1-T8	X-LARGE	12.0	32.0	76358330192	40	4	10	1.18	INDONESIA THAILAND MALAYSIA



100% PFAS FREE

\* Glove weights may vary (+/-0.3 grams)